

2D F-SENSOR

2D LAB
TECHPAP SAS ()



2D LAB

2D LAB

2D LAB

Techpap



Windows.

Samples	PM 1	Angle 1	Step 1	PM 2	Angle 2	Step 2	PM 3	Δ
1	4.20	114.14	0.76	0.58	101.86	1.05	0.16	8
2	4.20	114.14	0.76	0.58	101.86	1.05	0.16	8
3	4.20	114.14	0.76	0.58	101.86	1.05	0.16	8
4	4.20	114.14	0.76	0.58	101.86	1.05	0.16	8
5	4.20	114.14	0.76	0.58	101.86	1.05	0.16	8
6	4.20	114.14	0.76	0.58	101.86	1.05	0.16	8



TECHPAP SAS

CS 90251 - 38044 Grenoble cedex 9

Tél. +33 (0)4 76 51 74 75 - Fax +33 (0)4 76 42 05 04

www.techpap.com - techinfo@techpap.com